kasios® 

# JECTOS+®

**RADIOPAQUE Injectable**  
Synthetic Bone Substitute

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## Presentation

JectOS+ is an **injectable radiopaque synthetic bone substitute** made of Calcium Phosphate and radiopacifier (Zirconia-ZrO<sub>2</sub>).

The product is **easy** to mix and **to inject percutaneously**.

JectOS+'s excellent biocompatibility permits its osseointegration inside the vertebral body.

JectOS+ is sold packaged in a box containing:

- One vial of powder composed of the Calcium Phosphate and the **radiopacifier (ZrO<sub>2</sub>)**.
- One vial of liquid.
- Mixing accessories.
- A syringe to inject the product.
- A cannula.

## Indications\*

JectOS+ is an injectable bone substitute designed for **injection inside the vertebral body following a kyphoplasty on young patients with a stable and recent traumatic vertebral body fracture** (MAGERL classification, types A1 to A2).

**\* Caution: patients should wear a brace for the three days following surgery.**

## Contraindications

JectOS+ is contraindicated in unstable traumatic fractures, osteoporotic vertebral compression fractures and vertebral osteonecrosis.

## Chemical composition

The final product, post-chemical reaction, is a solid compound made of:

- TCP (Tricalcium Phosphate)
- DCPD (Dicalcium Phosphate Dihydrated)
- ZrO<sub>2</sub> (Zirconia)

The reaction is moderately exothermic, producing a transient rise in temperature (less than 45°C).

## Physical properties

Porosity: 40%

Pore size: < 5 microns

Setting time: 9 to 11 minutes (including **working time of 4 to 5 minutes\***)

Reaction temperature: 43°C to 45°C

**Compressive strength:** approximately **35 MPa**

*\* To increase working time, store JectOS+ in a refrigerator at +4°C.*

## Features and benefits

- **JectOS+ is highly radiopaque, so that progression in the vertebral body can be monitored during injection to avoid any risk of leakage.**
- JectOS+ **chemical reaction** is only **moderately exothermic**, thus minimizing shock reactions during injection, as do occur with PMMA type cements.
- The **35 MPa compressive strength** is higher than in cancellous bone and lower than in PMMA cement.
- **Easy injection** with a syringe or with an injection system.
- The **working time in excess of 4 minutes** is sufficient to fill completely the space left by the balloon.
- **JectOS+ stabilizes the vertebral body fracture, thereby alleviating or even suppressing fracture-induced pain.**
- JectOS+ restores vertebral body height and corrects kyphosis.

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**1. Pour entire quantity of liquid into bowl.**

Caution! Empty all the liquid in the vial into the bowl.

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**2. Add entire quantity of powder contained in the other vial into the bowl.**

Caution! Empty all the powder in the vial into the bowl.

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**3. Combine powder and liquid and stir vigorously for 30 seconds until mixture becomes smooth and of even consistency.**

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**4. Aspirate the liquid through the cannula connected to the syringe.**

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**5. From time to time, test a drop of the product to check for hardening.**

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**6. When JectOS<sup>+</sup> reaches a tooth-paste-like consistency, begin to inject into the vertebral body by connecting to the trocar already in situ.**

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### CAUTION!

**Do not inject JectOS<sup>+</sup> until its consistency is like toothpaste.**

As the chemical reaction is highly dependent on ambient temperature, setting times vary. For this reason we strongly recommend checking product hardening from time to time between minute 1 and minute 5.

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
Afshin Gangi, M.D., Professor  
 in cooperation with the Orthopædic Department  
 headed by Prof. J.-P. Steib  
 (Strasbourg - France)

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
- ◇ Mr B.L. - 27-year-old
- Injured while skiing (20.02.2006).
- Plain radiographs and CT scan showed A1 compression fracture of L1.
- No neurological deficit.
- Conservative treatment.
- Gradual worsening of symptoms.
- Referred to our department for kyphoplasty.
- Kyphoplasty done on 01.03.2006 (D9 of injury).
- The brace is removed three days after surgery.

**Imagerie fluoroscopique pendant la kyphoplastie**  
**Fluoroscopic images during kyphoplasty procedure**

1. Radiographie de contrôle.
  2. Introduction du trocart.
  3. Insertion du ballonnet de kyphoplastie.
  4. Gonflage du ballonnet.
  5. Ablation du ballonnet.
  6. Injection du JectOS+.
- 
1. Control radiograph.
  2. Trocar introduction.
  3. Insertion of kyphoplasty balloon.
  4. Balloon inflation.
  5. Balloon removal.
  6. JectOS+ injection.

**Comparaison entre imagerie TDM avant et après kyphoplastie**  
**Comparison between pre- and post kyphoplasty CT images**

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Sagittal reconstructed CT images before and after kyphoplasty with patient in prone position showing reduction of kyphotic deformity and improvement in vertebral height.

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
- ◇ Mrs J.A.V. - 31-year-old
- Injured during winter sport activity in Finland on 18.02.2006.
- Clinical presentation: backache with paresthesia.
- Imaging: burst fracture of L1 with retropulsed free fracture fragment in spinal canal.
- Referred back to France on 20.02.2006 for further management.
- Treatment plan: combined surgical-radiological intervention of open reduction internal fixation (T10-L2) with kyphoplasty of L1.
- Combined procedures performed on 23.02.2006 (D6 of injury).

**Radiographie standard et imagerie TDM (jour de l'accident)**  
**Plain radiographic and CT imagery (day of injury)**

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
**Fluoroscopic imagery during kyphoplasty**

Combined procedure\* on 6th day following injury.

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
1. Bilateral approach for kyphoplasty.
- 2 & 3. Balloon inflation.
4. JectOS+ injection into fractured L1.

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\* Type A3 fracture with retropulsed posterior free fracture fragment in spinal canal, surgically corrected, followed by kyphoplasty and filling of vertebral body with JectOS+ to avoid an anterior approach.

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Stabilization and healing of fractured L1 vertebra.  
 Restoration of vertebral height.  
 Absence of kyphotic deformity.

**Packaging**

Réf.

J+ 441Z20      10cc      Boîte de 1 - Stérile

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- 1 flacon de poudre
- 1 flacon de liquide
- 1 bol
- 1 spatule
- 1 seringue
- 1 canule

**Packaging**

Ref.

J+ 441Z20      10cc      Box of one - Sterile

Each box contains:

- 1 vial of powder
- 1 vial of liquid
- 1 mixing bowl
- 1 spatula
- 1 syringe
- 1 cannula

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